

ENGINEERING CHEMISTRY LABORATORY EXPERMENTS

B.Tech I YEAR

EXP...No. 1

TITRIMETRY


AIM: Estimation of hardness of water by EDTA method

APPARATUS: - Burette, Pipette, beakers, Standard flask, Conical flasketc.

CHEMICALS: - 0.1 M EDTA, Eriochrome, Blacke-T, water sample buffer solution.

THOERY:- Estimation of hardness by EDTA method is based on the following principles.

1. First, the indicator Eriochrome Black-T, which is a blue colored dye, forms an unstable complex with calcium or Magnesium ions in hard water at P^H of 9 to 10. The complex is wine red in color.



Unstable Complex (Wine Red)

2. As this solution is titrated against EDTA, the free Ca^{2+} and Mg^{2+} tones in water from stable metal ion EDTA $\rightarrow [\text{Mg/Ca-EDTA}]$ Stable Complex. Colorless

Once the free metal ions are complexed, the EDTA replaces Ca and Mg ions from the unstable indicator complex also, to form a stable complex, with the result, the indicator is set free. Since the free indicator is blue in colour at the above mentioned PH, the end point is the appearance of blue color



The temporary hardness is removed by boiling and after the removal of precipitate by filtration; the permanent hardness in the filtrate is determined by titration with EDTA as above.

Procedure:-

1. Rinse and fill the burette with EDTA solution.
2. Pipette out 25 ml of hard water into volumetric flask and add 20 ml of distilled water and 10 ml of buffer solution
3. Add 2-3 drops of Eriochrome Black-T indicator.
4. Titrate against 0.01M EDTA solution till the color of the solution changes from wine red to Blue.
5. Note the final reading and repeat to get three concordant readings.

Observation Table:

Sl.No	Volume of the solution taken in the titration flask (nl)	Burette Readings		Volume of the titrant used $V_2 - V_1$ nl
		Initial Reading (V_1 nl)	Final Reading (V_2 nl)	

→ Volume of hard water taken for each titration = 25 ml

→ Volume of 0.01 M EDTA used [from titration] = V_1 ml

Calculation: - 100 ml of 1 M EDTA = 100 gm of CaCO_3 = 1 mg of CaCO_3

1 ml of 0.01M EDTA = $100 \times \frac{1}{1000} \times \frac{1}{100\text{gm}}$ of CaCO_3 = 1 mg of CaCO_3

V ml of 0.01M EDTA = $I \times V$ Mg of CaCO_3 = x mg of CaCO_3

25 ml of hard water contains = x mg of CaCO_3

1000 ml of hard water contains = $x/25 \times 100$ mg of CaCO_3 = y Mg/L

RESULT: The Total hardness of water = y Mg/L

Precautions:-

1. The burette and pipette, conical Flask should be rinsed with distilled water.
2. The color change near the end point is very slow and this should be observed carefully.

EXPERIMENT NO.2

DETERMINATION OF Cu^{2+} THROUGH IODOMETRIC TITRATION

AIM: Determination of Cu percentage in Copper one (Brass)

APPARATUS: Burette, Pipette, Beakers, Titration flask

CHEMICALS: Copper One, Hypo, KI, Br_2 Solution, Acetic Acid, Nitric and HCl, H_2SO_4

INDICATOR: Starch solution

END POINT Blue to Colourless

PROCEDURE:

Preparation of Solution of ore

1. Take about 1 g of the finely ground copper ore in a beaker and warm with 15-20 ml of conc. HNO_3 until the copper passes into solution.
2. Evaporate to a volume of 5 ml and add 10 ml of conc. HCl and equal volume of 1: 1 H_2SO_4 .
Again evaporate till dense fumes of SO_3 appear indicating complete removal of nitrogen oxides.
3. Heat the mixture with 15 ml of bromine solution and an equal volume of water to complete oxidation of any As and Sb to the Pentavalent state.
4. Boil off excess of bromine and dilute the resulting solution to 100 ml in a measuring flask.

Titration:

1. Rinse and fill the burette with 0.1 N $\text{Na}_2\text{S}_2\text{O}_4$ solution
2. Measure 20 ml of the ore solution into a conical flask
3. Neutralize any free acid present by adding Na_2CO_3 solution drop by drop till a faint permanent precipitate remains even on shaking.
4. Add dil. Acetic acid drop wise until precipitate just dissolves
5. Add 1 gm of solid KI, cover the mouth of conical flask by filter paper and allow the mixture to stand for 2 to 5 minutes in the dark the solution becomes brown due to liberated iodine
6. Titrate the liberated iodine with hypo solution added from the burette, the brown colour of iodine becomes fainter as the addition by hypo solution proceeds and when only a faint yellow colour remains add 2 ml starch solution. This immediately form a deep blue iodo-starch complex
7. Now add hypo solution further, drop by drop when only milky white colour remains.

Observation Table:

Sl.No	Volume of the solution taken in the titration flask(ml)	Burette Readings		Volume of the titrant used (final-Initial reading((ml)
		Initial Reading	Final Reading	
1.				
2.				
3.				

Calculations:

(Ore Soln.) ($\text{Na}_2\text{S}_2\text{O}_3$)

$$N_1V_1 = N_2V_2$$

$$N_1 \times 20 = 0.1 \times V_2$$

$$N_1 = \frac{0.1 \times V_2}{20} \times \frac{63.5}{\text{Strength of Solution}} \times 100$$

$$N(\text{Cu}^{2+}) = \frac{V_2}{200} \times \text{Strength of } \text{Cu}^{2+} = \frac{V_2}{200} \times 63.5 \text{ g/L}$$

$$\text{Percentage of Cu in ore solution} = \frac{V_2}{200} \times \frac{63.5}{\text{Strength of Solution}} \times 100$$

RESULT: The amount of Cu present in copper ore solution is

EXPERIMENT NO.4

COLORIMETRY

ESTIMATION OF COPPER BY COLORIMETRIC METHOD

AIM: Estimation of Copper by Colorimetry

APPARATUS: - Colorimeter, Test tubes, Burettes.

PRINCIPLE: - Colorimeter measures the optical density of an absorbing substance where optical

density (O.D) is defined as $O.D = \log \frac{I_o}{I} \rightarrow 1$

Where I_o = Intensity of incident light

I = Intensity of transmitted light

As per beers law, optical density of an absorbing substance is related to the concentration by the equation.

$$O.D = E.C.l$$

$$O.D = (E.l).C \rightarrow 2$$

Where 'C' is the concentration of the substance, l is the path length, which represents the width of the cell used and is constant for a given cell used, E is the molar absorption coefficient and is a constant for given substance. Equation 2 may be written as $O.D. \propto C \rightarrow 3$

Equation 3 represents the quantitative form of Beer's law . if the optical density of a substance is determined at varying concentration. A plot of O.D.vs C gives a straight line.

PROCEDURE:- take the sample solution of CuSo₄ and prepare the following 10 sample solutions in test tubes as 1 to 10

1. 1 ml 0.1 m CuSo₄ + 9 ml Distilled water
2. 2 ml 0.1 m CuSo₄ + 8 ml Distilled water
3. 3 ml 0.1 m CuSo₄ + 7 ml Distilled water
4. 4 ml 0.1 m CuSo₄ + 6 ml Distilled water
5. 5 ml 0.1 m CuSo₄ + 5 ml Distilled water
6. 6 ml 0.1 m CuSo₄ + 4 ml Distilled water
7. 7 ml 0.1 m CuSo₄ + 3 ml Distilled water
8. 8 ml 0.1 m CuSo₄ + 2 ml Distilled water
9. 9ml 0.1 m CuSo₄ + 1 ml Distilled water
10. 10 ml 0.1 m CuSo₄ + 0 ml Distilled water

The ten sample solutions prepared above have a varying concern from 0.01m to 0.1m, choose the filter in the colorimeter with maximum absorbance. Tabulate the result of filter and O.D with a given CuSO_4 sample solutions.

After selecting filter, determine the O.D. of the above mentioned ten sample solutions and tabulate the results.

Observation:

Table-1

Filter No	λ -Range (nm)	Peak (nm)	Optical Density(OD)
1			
2			
3			
4			
5			
6			
7			
8			

Table-II

Selected Filter-----

Sl.No	Volume of CuSO_4 (ml)	Volume of H_2O (ml)	Concentration of CuSO_4	Optical Density
1	10	0		
2	9	1		
3	8	2		
4	7	3		
5	6	4		
6	5	5		
7	4	6		
8	3	7		
9	2	8		

10	1	9		
----	---	---	--	--

GRAPH: - Plot a graph of O.D Vs [CuSo₄]. A straight line passing through origin is obtained the slope of which gives El, E is the molar absorption coefficient.

$$\text{Slope} = \frac{y_2 - y_1}{x_2 - x_1}$$

$$\text{slope} = El \Rightarrow E = \frac{\text{slope}}{l}$$

$$E = -\text{mole}^{-1}\text{cm}$$

RESULT: - Amount of copper present in given CuSo₄ solution is _____grams?

EXP.NO.5

CONDUCTOMETRY

Aim: - To determine strength of given HCl Solution Conductometrically

APPARATUS: - Conductivity meter, conductivity cells, Beakers, Pipette, Burette

CHEMICALS:- Oxalic acid Solution (0.1N), NaOH Solution, HCl, Phenolphthalein Indicator.

Preparation of Standard Oxalic acid Solution (0.1N)

About 1.6g of Oxalic acid dehydrate is weight accurately and dissolved in small amount of distilled water and made up to the mark in 250 ml volumetric flask.

PRINCIPLE: - The conductivity of HCl is very high due to mobility of H^+ ions in it. When a solution of strong alkali is added to the solution of strong acid, the highly conducting H^+ ions are replaced by Na^+ ions is considerably smaller than H^+ and hence addition of alkali to the acid always accompanied by a decrease in conductivity. When all the H^+ are replaced that is after the end point, further addition of NaOH results in increase in conductivity. From the plot of the conductive Vs Volt NaOH added, the end point can be calculated

PROCEDURE: - 20 ML OF Oxalic acid is pipette out into a conical flask and is diluted to 30 ml with distilled water. To the solution, two drops of phenolphthalein indicator is added and the titrated against sodium hydroxide taken in a burette. The endpoint colorless to pale pink. The experiment is repeated till correct reading are obtained

Determination of Strength of Given HCl Solution:-

The given unknown HCl solution is 100 ml volumetric flask is made up to the mark and then 10ml of given unknown HCl Solution is taken in a 100 ml beaker. To this 40ml of distilled water is added, the electrode of the cell is dipped in the solution. The conductivity of solution is noted before adding the alkali. The standardized sodium hydroxide is added from a burette with 1ml of increment and by shaking thoroughly, the contents of the beaker the conductivity values are noted. The values of observed conductivities are plotted against volume of sodium hydroxide. The concentration of hydrochloric acid is determined using the end point from the graph.

PRECAUTIONS: - after 2 each addition of titrant from the burette the solution should be thoroughly stirred and then the reading should be taken.

RESULT: - The Strength of Hydrochloric Acid is _____ g/lit.

Calculation:-

Standardization of Sodium Hydroxide Solution

S.No	Volume of Oxalic Acid (V ₁ ml)	Burette Readings		Volume of NaOH Solution (V ₂ =F-I)
		Initial	Final	
1	20			
2	20			
3	20			

WE know that $V_1N_1 = V_2N_2$

V₁ is Volume of Oxalic Acid = 20ml

N₁ is Normality of Oxalic Acid= 0.1N

V₂ is Volume of Sodium Hydroxide= - ml

N₂ is Normality of Sodium Hydroxide – N

$$N_2 = \frac{V_1 N_1}{V_2}$$

∴ Normality NaOH = -----N

Determination of Strength of HCl:-

Sl.No	Volume of NaOH	Conductivity (ohm ⁻¹)
1		
2		
3		
20		
40		

Volume of NaOH Solution from graph Ph V₂ = -ml

GRAPH:

Normality of NaOH Solution N₂= -N

Volume of Hydrochloric acid solution V₃= 10ml

Normality of HCl in mixture N₃= -----N.

$$V_2 N_2 = V_3 N_3$$

$$N_3 = \frac{V_2 N_2}{V_3}$$

Strength of HCl Solution = Normality of HCl Solution \times eg. Wt of HCl

$$= N \times 36.5$$

$$= n \text{ _____ g/Lit}$$

EXP.NO.6

POTENTIOMETRY

TITRATION OF STRONG ACID Vs STRONG BASE

AIM:- To Determine the concentration of a strong acid by a Potentiometry titration given a strong base of known concentration.

APPARATUS: - Potentiometer, Pt-electrode, Calomel electrode, Salt bridge, 150 ml beakers, Stirrer, Pipette, Burette.

CHEMICALS: - 0.1m NaOH, 0.1m HCL, Quinhydrone solid, KCl Solution.

PRINCIPLE: - The quinhydrone electrode is prepared by taking saturated solution of quinhydrone in a known quantity of strong acid whose concentration is to be determined. This electrode is combined with the calomel electrode to make the cell. Pt, Hg, $\text{Hg}_2\text{Cl}_2(\text{s})/\text{Saturated}/\text{H}^+/\text{QH}_2/\text{Pt}$

In acid medium quinhydrone exist in two forms.

As NaOH is added to a solution of H^+ ions with quinhydrone the H^+ ions are consumed steadily and the relative concentration of QH_2 and Q changes with the titration. At the end point where all the ions are consumed, QH_2 gets completely converted to Q. the electrode potential of the electrodes is given as per the nearest equation.

$$E_{\text{QH}_2 / \text{Q}} = E^0_{\text{QH}_2 / \text{Q}} - \frac{0.059}{n} \log \frac{a_{\text{Q}} \cdot a_{\text{H}^+}}{a_{\text{QH}_2}}$$

$$= E^0_{\text{QH}_2 / \text{Q}} - \frac{0.118}{n} \log a_{\text{H}^+} \quad [\because a_{\text{Q}} \text{ and } a_{\text{CH}_2} = 1]$$

At the endpoint, which $a_{\text{H}^+} \rightarrow 0$ there is a sudden fall in the cell emf $E_{\text{cell}} = E_{\text{calomel}} + E_{a_{\text{H}_2}/\text{Q}}$

Since calomel is constant the variation of E_{cell} is a measure of the variation of $E_{aH_2/Q}$. The endpoint is the titration is the point where there is a maximum fall in potential with a small addition of NaOH.

Procedure: - Pipettes out 10 ml of the given HCl solution into a clean 100ml beaker add a pinch of quinhydrone and stir it well to make a saturated solution of quinhydrone is H^+ . dip a platinum electrode. Take a second beaker with saturated KCl solution in it. Dip the calomel electrode in it. Connect two beakers with a salt bridge; connect the terminals of the potentiometer with the two electrodes. Note the emf of the cell. fill a burette with 0.1 NaOH, add 0.5ml of NaOH into the beakers containing acid, stir well for equilibrium to establish. Note the emf of the cell once again. The variation of emf by successive addition of 0.5 ml of NaOH. Each time will be about 10-20mv, as the endpoint approaches, the fall in the emf would increase from 10-20 mv to 40-50 mv. At this stage add NaOH in batches of 0.2 ml and note emf after each addition. At the end point the emf falls suddenly by about 100-150 mv. Continue the titration with 0.2 ml of NaOH and note 5-6 emf values. Now the variation in emf gets smaller and takes 4-5 emf readings, after the endpoint by adding batches of 1ml of NaOH each time, tabulate the results.

GRAPH: - Plot a graph of emf Vs volume of NaOH added and a second graph of $\Delta E / \Delta V$ Vs Volume of NaOH added. Note the point of Neutralization from both the graph and calculate the conc of given HCl using the conc of given NaOH as 0.1M

RESULT: - At the Neutralization point 10 ml of given HCl =? And conc of given HCl=_____

Sl.No	Volume of Alkali (ml)	Emf (mv)	$\xrightarrow{\text{Poly'n}} \frac{\Delta E}{\Delta V}$
1			
2			
3			

. . . 30			
-----------------------	--	--	--

CALCULATION: -

HCl

NaOH

$$V_1 = 10ml$$

$$V_2 = \text{_____ (From Graph)}$$

$$N_1 = ?$$

$$N_2 = 0.1$$

$$V_1 N_1 = V_2 N_2 \Rightarrow N_1 = \frac{V_2 N_2}{V_1}$$

∴ Concentration of given HCl (N_1) = ____N.

EXPERIMENT – 9

IDENTIFICATION OF FUNCTIONAL GROUPS IN ORGANIC COMPOUNDS

SYSTEMATIC PROCEDURE

EXPERIMENT

OBSERVATION

INFERENCE

PHYSICAL CHARACTERISTICS:

a) Physical state

i) Liquid (or) Solid)

ii) Amorphous (or)

crystalline

b) Odour

i) Pale yellow

Nitro compounds

ii) Reddish orange

Nitro amines

	iii) Brownish yellow	Aniline
	iv) Pale pink	Phenol
	v) Colourless	Carbohydrates, aldehydes, ketones, carboxylic acids, esters and some simple hydrocarbons.
c) Odour	i) Odourless	Carbohydrates, aromatic acids
	ii) Bitteralmonds	Nitrobenzene, benzaldehyde are absent
		Aromatic amines
	iii) Fishy smell	Acetamide
	iv) Mouse like smell	Formaldehyde, formic acid
	v) Pungent smell	Esters
	vi) Pleasant fruites	Phenols
	(Phenolic) Carbolic smell	Acetic acid
	Vinegar like Cinamon like	Cinnamaldehyde

d) Melting point/ Boiling point

II. Solubility : Take a small amount of compound in a test tube and add 3-ml of solvent

a) Water	Soluble in cold water	i) Carbohydrates, glucose, Fructose
		ii) Amides, urea acetamide
b) 5% NaOH	Soluble in 5% NaOH and regerated with HCl	May be week acid like phenols

c) 5% NaHCO ₃	Soluble in 5% NaHCO ₃ and regenerated with HCl	May be strong acids, benzolic, salicylic acids.
d) 5% HCl	Soluble in 5% HCl regenerated with NaOH.	May be amines like aniline
e) Ether	Soluble in ether	May be neutral compounds like, chlorobenzene, Aldehydes, Ketones, Nitrobenzene etc.

III. Ignition Test :

Heat a small amount Of given compounds On spatula	1) Burns with non-luminous and non-smoky flame.	Aliphatic compound
	2) Burns with Luminous And smokes flame	Aromatic compound
	3) Burns with sugar smell	Carbohydrates

DETECTION OF FUNCTIONAL GROUPS

TEST FOR PHENOLS

i) With neutral FeCl₃ test :

The compound to dissolved	Violet (or) Blue	Phenol is present
In alcohol and add 1 (or) 2 Drops of 5% FeCl ₃	colour is observed	



ii) Liebermann's nitroso reaction:

Take the compound in a dry Test tube add sodium nitrite And conc sulphuric acid(1ml) Mix well and heat gently	A blue colour is obtained the solution turns red on dilution with water and blue on basification with dil Sol of NaOH	Phenol is present
---	---	-------------------

1.NaHCO₃ Test

$$R-COOH + NaHCO_3 \rightarrow RCOONa + H_2O + CO_2 \uparrow$$
$$R-COOH + R^1OH \xrightarrow[Heat]{Conc. H_2SO_4} RCOOR^1 + H_2O$$

Test for Carbonyl Compounds:

Dissolve the carbonyl compound	a yellow or orange colour	may be a
(100 mg (or) 1-2 drops) in ethanol (2-3ml)	ppt is formed	carbonyl compound
To this add an alcoholic solution of		
2,4 – di nitro-phenyl hydrazine (2 ml)		
and Shake the mixture well.		

ii) Tollen's Test:


To the Tollen's reagent (1ml)
Add a solution of the aldehyde
(2 drops (0.1) 50 mg) dissolved
in aldehyde free alcohol
(2-3 ml) and warm the solution
in a hot water bath.

a grey black ppt (or) silver
mirror deposits on the inner
walls of a test tube May be aldehyde is
present

iii) Fehling Test:

Add Fehling solution (2-3 ml)
By mixing equal amounts of
Fehling A + Fehling B to the
Organic compound and warm the mixture

Red precipitate of
cuprous oxide is formed may be aldehyde is
present



Fehling Solution

Red

iv) Schiff's Test

To a solution organic
Compound add Schiff's
Reagent (2-3 drops) and
Shake the contents

Violet (or) purple colour (pink)
is observed


may be aldehyde is
present

v) Iodo form test:

Dissolve the given compound
(2-3 ml (or) 100 mg) in water
(2-3 ml) in a test tube and add
2-3 ml NaOH (10%). To this add
A saturated solution of iodine in KI
with stirring until the dark colour of
iodine persists. Heat the solution in a
boiling water bath for 1-2 minutes.

A yellow ppt Iodo form
is observed

may be ketone



Test for Carbohydrates:

i). Molisch's Test:

Add an alcoholic solution of
 α -Naphthol (10%) to an aqueous
Solution of substance. Then add
Conc H_2SO_4 (1ml) carefully along
the sides of the test tube allow to stand
for 2 minutes.

A violet ring is observed at the
junction of two layers

may be carbohydrate

ii).Barfoed Test:

Heat the aqueous solution of the compound (1ml) with barfoed reagent (1ml) in boiling water bath for 1-2 minutes

i) red colour precipitate of Cu ₂ O is formed	carbohydrate is present present
--	---------------------------------

iii).Test with Benedict solution:

Add Benedict solution (2ml) to the Aqueous solution of organic substance (reducing sugar)


Red ppt is observed	Carbohydrate is present
---------------------	-------------------------

IV). TEST FOR [1⁰, 2⁰, 3⁰] AMINES:

1. Carbylamine test:

To the compound add few ml Of chloro form and then add 2 ml of alcoholic KOH; mix well and gently warmed.

a foul smell is observed	may be 1 ⁰ amine
--------------------------	-----------------------------



2. Azodye Test:

Take 2 ml of compound, add 2 ml Of HCl, coll in ice, and then add 2 ml of ice cold 10% Aq NaNO₂ solution. Add to it cold solution of 0.4 gms of 2-Naphthol in 4 ml of 5% NaOH solution is added.

orange red dye is formed	presence of 1 ⁰ amine
--------------------------	----------------------------------

EXPERMENT NO.10

Preparation of Aspirin

AIM : To Prepare Aspirin from salicylic Acid

APPARATUS: Conical flask, Water bath, beakers, glass rod

CHEMICALS: Salicylic acid, acetic anhydride, sulphuric acid, alcohol

PROCEDURE:

1. Take 5 gm of salicylic acid in a 100 ml conical flask. To this add about 10 ml of acetic-anhydride and 1-2 ml of con.H₂SO₄.
2. 2Shake the contents to 800C, because of exothermic reaction. Maintain the tem at 60-70⁰ for about 15 minutes, keeping on a water bath.

3. Allow the solution to come to room temperature and then pour it in 100 ml cold water taken in a 500 ml beaker with stirring (water is added to destroy the excess acetic anhydride which gets converted to acetic acid).
4. To induce crystallization scratches the sides of the flask with glass rod.
5. Filter the solid and wash it with cold water. Press it in between the folds of filter paper.
6. Spread the solid on the filter paper to dry
7. To obtain colorless crystals, recrystallization can be done using equal volumes of ethanol and water solvent system.

RESULT: - Yield=

Precautions:-1 Take acetic anhydride in excess as it acts as acetylating agent as well as Solvent.

2. Make sure that all the salicylic acid is dissolved.

CHEMICAL EQUATION:

Experiment No.11

KINETICS:

To determine the rate constant of hydrolysis of methyl acetate catalysed by an acid also the activation energy.

AIM: To determine the rate constant of hydrolysis of methyl acetate catalysed by two different concentrations of acids.

APPARATUS: - Conical Flask, Beakers, Reagent bottles, Burette, Pipette.

CHEMICALS: - Methyl acetate, Hydrochloric acid (1M, 2M)


Oxalic acid solution (0.2N), NaOH, Phenol phthalein indicator

Preparation of Standardized oxalic acid solution (0.2N)

About 3.15 g of oxalic acid di hydrate is weighed accurately and dissolved in small amount of distilled water and made up to the mark in 250ml of volumetric flask.

$$N = \frac{W}{E \cdot V} \times \frac{1000}{250}$$

Principle: - The hydrolysis of methyl acetate in presence of an acid may be represented as



The rate of reaction is given by

$$\frac{dx}{dt} = \frac{2.303}{t} \log\left(\frac{a}{a-x}\right) \rightarrow (\text{First Order})$$

As acetic acid is produced during the reaction, the reaction may be studied by the titration of unknown concentration of reaction mixture with a standard alkali at suitable time.

Standardization of NaOH Solution-

20 ml of Oxalic acid is pipette out into a conical flask and is diluted to 30 ml with distilled water. To the solution two drops of phenolphthalein indicator is added and then titrated against sodium hydroxide taken in a burette. The end point is colorless to pale pink. The experiment is repeated till concurrent readings are obtained.

PROCEDURE:

From standardized 1N and 2N HCl, 10 ml of 1N HCl taken in a iodination flask and 10 ml of pure methyl acetate is taken in a test tube. Both the iodination flask and test tubes are placed in water bath to bring them into room temperature. Measure exactly 5 ml of methyl acetate and add to 90 ml 1N HCl in a iodination flask. When exactly half of the amount in the pipette is been discharged, the stop watch is started. Mixture is shaken well and maintained at a constant temperature. Now titrated with NaOH solution + using phenolphthalein indicator. The volume of alkali required for titration is noted and is taken as " V_0 ". Similarly 10 ml of reaction mixture is taken out in successive intervals of time 10, 20, 30, 40, 50 minutes. To obtain complete hydrolysis the reaction mixture kept in a water bath maintained at for about half an hour.

After the hydrolysis completed and it is cooled at room temperature and titrated the 10 ml of reaction mixture with NaOH solution and is taken as " V_∞ ".

Method of Calculation:-

The initial volume of titrate; NaOH is V_0 at T_0 °C is proportional to the amount of HCl present in the reaction mixture at Zero time. When there is an acetic acid present as the hydrolysis proceeds $V_\infty - V_0$ directly proportional to initial concentration of ester (a), at successive intervals of time $V_\infty - V_t$ is proportional the concentration of esters at time 't' (a-x).

$$K_1 = \frac{2.803}{t} 10g \left(\frac{V_\infty - V_0}{V_\infty - V_t} \right) 80C^{-1}$$

Same procedure and calculation is repeated for 2N HCl solution.

Relative strength of two acids is given by $\frac{K_2}{K_1} =$

Calculation of Activation Energy:-

$$\log \frac{K_2}{K_1} = \frac{Ea}{2.303 \times R} \times \frac{T_2 - T_1}{T_2}$$

$$Ea = 10g \frac{K_2}{K_1} \times 2.3.3 \times R \times \frac{T_2}{T_2 - T_1}$$

Result:-Rate Constant for hydrolysis of methyl acetate by 1N HCl are $K_2 = \text{-----}$, $K_2 = \text{----}$

Activation Energy $Ea = \text{-----}$

Observation Table:-

Table-1 Standardization of NaOH:-

Sl.No	Vol of Oxalic acid (V1M1)	Burette Readings		Vol of NaOH (V2) V2= F-I
		Initial	Final	
1	20			
2	20			
3	20			

Vol of Oxalic acid Solution $V_1 = 20$ ml

Normality of Oxalic Acid $N_1 = 0.2$ N

Vol of NaOH Solution $V_2 = \text{----}$ M1

Normality of NaOH $N_2 = \text{-----}$ N

$$N_1 V_1 = N_2 V_2 \Rightarrow N_2 = \frac{N_1 V_1}{V_2} \rightarrow$$

\therefore Normality of NaOH $N_2 = \text{-----}$ N

Calculations for determine rate constants

For 1N Hcl

Sl.No	Time (t) Minutes	Vol of NaOH (Vt)	$V_{\alpha} - V_t$	$\frac{V_{\alpha} - V_0}{V_{\alpha} - V_t}$	$10g \frac{V_{\alpha} - V_0}{V_{\alpha} - V_t}$	$K_1 = \frac{2.303}{t} \log \frac{V_{\alpha} - V_0}{V_{\alpha} - V_t}$
1	0					
2	10					
3	20					
4	30					
5	40					
6	50					
7	60					
8	V_{α}					
9						

Average $K_1 =$

$$V_0 = \text{---} ml$$

$$V_{\alpha} = \text{---}$$

$$V_{\alpha} - V_0 = \text{---} ml$$

From Graph

$$Slope = \frac{y_2 - y_1}{x_2 - x_1} = \text{---}$$

$$Slope = \frac{K_1}{2.303}$$

$$\Rightarrow K_1 = 2.303 \times Slope$$

GRAPH:

For 2N HCl:

Sl.No	Time (t) Minutes	Vol of NaOH	$V_{\alpha} - V_t$	$\frac{V_{\alpha} - V_0}{V_{\alpha} - V_t}$	$10g \frac{V_{\alpha} - V_0}{V_{\alpha} - V_t}$	$K_2 = \frac{2.303}{t} \log \frac{V_{\alpha} - V_0}{V_{\alpha} - V_t}$
-------	---------------------	----------------	--------------------	---	---	--

		(Vt)				
1	0					
2	10					
3	20					
4	30					
5	40					
6	50					
7	60					
8	V_{α}					
9						

Average $K_2 =$

$$V_0 = \dots ml$$

GRAPH:

$$V_{\alpha} = \dots$$

$$V_{\alpha} - V_0 = \dots ml$$

$$Slope = \frac{y_2 - y_1}{x_2 - x_1} = \dots$$

$$\text{From Graph } Slope = \frac{K_2}{2.303}$$

$$\Rightarrow K_2 = 2.303 \times Slope$$

EXPERIMENT -12

AIM: - Preparation of Thiokol Rubber

Operators: - Beakers, Glass Rod, Funnel, Filter Paper

Chemicals: - Ethylene Dichloride, Na_2S_2 (Sodium Poly Sulphide)


Procedure: -

1. Place about 10 ml. of Ethylene dichloride solution in 100 ml. beaker
2. Add about 5 gr. of sodium polysulphide and Heated to form poly sulphide rubber (solid)
3. Wash the above solid with water and dry it in the folds of filter paper
4. Calculate the yield of the product

Result:

The yield of Thiokol Rubber is _____ gr.

CHEMICAL EQUATION:



(Ethylene Chloride) (Sodium Poly Sulphide) (Poly Sulphide Rubber (or) Thiokol)